Rapid Exchange PTCA Balloon Dilatation Catheter

ACROSS HP™ **ACROSS CTO™_{RX}** **ACROSS CTO™_{ST}**

INSTRUCTIONS FOR USE	English	2
GEBRAUCHSANLEITUNG	Deutsch	4
MODE D'EMPLOI	Français	6
ISTRUZIONI PER L'USO	Italiano	8
INSTRUCCIONES DE USO	Español	10
MODO DE EMPREGO	Português	12
GEBRUIKSAANWIJZING	Nederlands	14
BRUGSVEJLEDNING	Dansk	16
BRUKSANVISNING	Norsk	18
BRUKSANVISNING	Svenska	20
INSTRUKCJA UŻYCIA	Polskie	22
NÁVOD K POUŽITÍ	Čeština	24
NÁVOD NA POUŽÍVANIE	Slovenský	26
HASZNÁLATI UTASÍTÁS	Magyar	28
INSTRUCȚIUNI DE UTILIZARE	Română	30
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	Ελληνικά	32
РУКОВОДСТВО ПО ПРИМЕНЕНИЮ	Русский	34
INSTRUKSI PENGGUNAAN	Indonesiai	36
GRAPHICAL SYMBOLS		38

Acrostak (Schweiz) AG
Stegackerstrasse 14
8409 Winterthur, Switzerland
Phone +41 52 233 95 51,
clinical@acrostak.com www.acrostak.com

QAdvis EAR AB
Ideon Science Park
Scheelevägen 17
SE-223 70 Lund, Sweden

CE 2797

1 DEVICE NAME

The device brand names are **ACROSS HP™ / ACROSS CTO™ RX / ACROSS CTO™ ST**. The generic device name is PTCA balloon dilatation catheter.

2 DEVICE DESCRIPTION

The PTCA balloon dilatation catheter is designed for dilating stenotic atherosclerotic lesions in coronary arteries or bypass grafts (PTCA). The catheter's dilating element is a balloon near the distal tip. The balloon has either one or two radiopaque marker band(s) (see package label for specific model). The catheter has two visual markers proximal on the shaft located approximately 90 cm (brachial access) and 100 cm (femoral access) from the distal tip that indicate, approximately, the exit of the catheter tip from the guiding catheter. The shaft has a distal port (hole) approximately 23 cm from the distal tip that accesses the central guidewire lumen. The guidewire lumen begins at the distal port and ends at the distal tip. The maximum guidewire diameter that may be used is printed on the package label. At the proximal end the catheters have a hub that is the balloon inflation port. The balloon is inflated by injecting diluted contrast medium through this hub. The proximal hub does not have a guidewire lumen. The guidewire exits the guidewire lumen at the distal port and travels alongside the catheter shaft and hub for the proximal portion of the catheter. This design allows insertion and removal of the catheter without extension of the guidewire. The compliance chart shows how the balloon diameter increases as inflation pressure increases; the chart also includes the pressure at which the balloon inflates to its nominal diameter (see package label for specific model). The use of a guiding catheter (dimension see package label for specific model) is recommended.

3 INDICATIONS

The PTCA catheters are indicated for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion.

4 CONTRAINDICATIONS

- Unprotected left main coronary artery.
- Coronary artery spasm in the absence of a significant stenosis.

5 WARNINGS

- This product is intended for single use only. Do NOT reprocess, resterilise and/or reuse. Processing and/or resterilisation measures can harm the characteristics of the catheter (for example, the inflation behavior, deflation behavior, leak tightness or tensile strength) due to the chemical, thermal and mechanical stresses. Furthermore, after reprocessing there is an increased risk of contamination of the product as well as infection or cross-infection of the patient.
- To reduce the potential for vessel damage, the inflated diameter of the balloon should not exceed the diameter of the vessel proximal and distal to the stenosis.
- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is

met during manipulation, determine the cause of the resistance before proceeding.

- Balloon pressure should not exceed the rated burst pressure. Use of a pressure monitoring device is recommended to prevent over pressurization.
- PTCA should only be performed at hospitals where emergency coronary artery bypass graft surgery can be immediately performed in the event of a potentially injurious or life-threatening complication.
- Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate the balloon.
- Do not use the catheter after the "Use by" date specified on the package.

6 PRECAUTIONS

- Ensure that the catheter is suitable for the intended procedure.
- The catheter should be used only by physicians trained in the performance of PTCA.
- Utilize appropriate anticoagulant/antiplatelet therapy and coronary vasodilator agents during procedure.
- During withdrawal of the PTCA catheter, hold a saline soaked gauze around the exposed catheter shaft and pull the catheter through the gauze to remove any excess contrast medium.
- Prior to insertion or withdrawal of the PTCA catheter, wipe the guidewire with saline soaked gauze to remove any contrast medium.
- Do not bend or kink the hypotube or the attached catheter's hub (balloon inflation port). Do not expose the catheter to organic solvents (e.g. alcohol).
- Do not autoclave. Exposure to temperatures above 54°C may damage the catheter.

7 POSSIBLE ADVERSE EFFECTS

Potential adverse events include, but are not limited to, the following:

- Dissection or perforation of the coronary artery
- Injury, rupture, intimal tear or other damage to the coronary artery
- Thrombosis, sudden and total occlusion of the coronary artery or bypass graft
- Unstable angina
- Arterial spasm
- Arteriovenous fistula
- Hypo/hypertension
- Hemorrhage or hematoma
- Allergic reactions
- Infections
- Acute myocardial infarction
- Arrhythmias
- Ventricular fibrillation
- Embolism
- Restenosis of the dilated vessel
- Death

8 HOW SUPPLIED

STERILE: Sterilized with ethylene oxide gas. Non-pyrogenic. Do not use if sterile package is opened or damaged.

CONTENTS:

- One (1) PTCA balloon dilatation catheter
- One (1) Instructions for use manual
- One (1) Compliance card

STORAGE: Store at room temperature, in a dry, dark place.

DISPOSAL INSTRUCTIONS: After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

9 OPERATOR'S INSTRUCTIONS

9.1 Inspection Prior to Use

Carefully inspect the sterile package before opening. Do not use if any defects or kinks are noted or after the "Use by" date. If the integrity of the sterile packaging has been compromised before exceeding the product's "Use by" date, contact your local representative for return information.

9.2 Materials required (not included in device package)

Appropriate guiding catheter(s) (see product label for dimensions)

- 20 ml (cc) syringes
- Heparinized normal saline solution (sterile)
- ≤ 0.014 inch (0.36 mm) guidewire
- Rotating hemostasis valve
- Diluted contrast medium 1:1 with sterile normal saline solution
- Inflation device
- Three-way stopcock
- Torque device
- Flushing needle

9.3 Preparation

9.3.1 Packaging Removal

1. Carefully remove the catheter from its protective dispenser.

9.3.2 Catheter Preparation

1. Attach a three-way stopcock to catheter's hub (balloon inflation port).
2. Attach a 20 ml syringe to the stopcock. Open the stopcock and induce negative pressure, pulling the syringe plunger as far back as possible without dislodging it from the syringe barrel.
3. Close the stopcock on the inflation port.
4. Remove the syringe and purge the air.

5. Reattach the syringe. Open the stopcock and again induce negative pressure, pulling the syringe plunger as far back as possible without dislodging it from the syringe barrel.

6. Close the stopcock on the inflation lumen and remove the syringe.

NOTE: To ensure air contained in the balloon and inflation lumen is removed, it is recommended that negative pressure be induced twice. Open the stopcock only when the syringe is in place with negative pressure induced.

7. Prepare an angioplasty inflation device with diluted contrast medium and attach it to the stopcock.

NOTE: Purge the angioplasty inflation system and the stopcock of air prior to opening the stopcock to the PTCA catheter.

CAUTION: Nonionic contrast medium has different viscosity and precipitation levels than does the ionic type, which may prolong inflation/deflation times.

8. With tip down, orient catheter vertically.
9. With the protective tubing and packaging stylet in place, apply positive pressure and open the stopcock to allow the contrast solution to flow into the inflation lumen and balloon. Verify the integrity of the catheter by pressurizing it to 5 atm (507 kPa).
10. Induce negative pressure and close the stopcock. Purge any air from the angioplasty inflation device through the stopcock.

NOTE: If air is present in the inflation device or PTCA catheter, repeat steps 9 and 10 until air has been removed.

11. Leave the catheter under negative pressure until it is ready to be used.
12. Remove the protective tubing with packaging stylet from the balloon.

9.3.3 Guidewire Lumen Flush

1. Flush balloon catheter guidewire lumen with heparinized saline using flushing needle (not included).

9.4 PTCA Procedure

1. Prepare the vascular access site according to standard PTCA practice.
2. Attach a hemostasis valve to the luer port of the guiding catheter positioned within the vasculature.
3. Insert the guidewire into the guiding catheter through the hemostasis valve. Under fluoroscopy, position the guidewire across the stenosis according to accepted PTCA techniques.
4. Insert the proximal end of the guidewire into the distal tip of the PTCA catheter. The guidewire will exit through the distal port (guidewire exit).

NOTE: Before inserting the PTCA catheter, wipe the guidewire with a saline soaked gauze to remove any excess contrast medium.

CAUTION: Care should be taken when inserting the PTCA catheter into the hemostasis valve in order to avoid kinking.

5. Advance the PTCA catheter to the distal end of the guiding catheter.

NOTE: The proximal exit markers on the catheter shaft may be used to estimate when the distal end of the PTCA catheter has reached the distal end of the guiding catheter.

CAUTION: Fully deflate the balloon by inducing negative pressure with the inflation device whenever advancing or withdrawing the PTCA catheter.

CAUTION: Do not advance or withdraw the PTCA catheter within the coronary vasculature unless it is preceded by a guidewire.

CAUTION: If a Tuohy-Borst type adjustable hemostasis valve is used, avoid overtightening since this action may restrict the flow of contrast medium in and out of the inflation lumen, thereby slowing inflation/deflation.

6. Continue under fluoroscopy and use the radiopaque marker band(s) to position the balloon within the stenosis.

7. Close the hemostasis valve and dilate the stenosis using an accepted coronary angioplasty technique.

NOTE: Do not exceed the rated burst pressure labeled.

NOTE: Balloon inflation should be performed with the guidewire extended beyond the catheter tip. It is strongly recommended that the guidewire, PTCA catheter, or both, remain across the stenosis until the procedure is complete.

9.5 Catheter exchange procedure

1. Ensure that the balloon is fully deflated.
2. Loosen the hemostasis valve.
3. Hold the guide wire and the hemostasis valve in one hand, while grasping the balloon catheter shaft with the other hand. Maintain the guide wire's position in the coronary artery by holding the guide wire stationary while carefully pulling the catheter out of the guiding catheter

CAUTION: Monitor the guidewire position under fluoroscopy during the exchange.

4. Pull the catheter until the guidewire exit is reached. Then proceed the same way as with an over the wire system until the catheter has been totally removed from the guiding catheter. Tighten the hemostasis valve.
5. Completely remove the PTCA catheter from the guidewire.
6. Prepare and insert the next catheter (according to the manufacturer's *instructions for use*).

9.6 Removal Procedure

1. Ensure balloon is fully deflated.
2. Fully open hemostasis valve.
3. While maintaining guidewire position and negative pressure on inflation device, withdraw the PTCA catheter.

NOTE: Before withdrawing the PTCA catheter, wipe the guidewire with a saline soaked gauze to remove any excess contrast medium.

NOTE: During withdrawal of the PTCA catheter, hold a saline soaked gauze around the exposed catheter shaft and pull the PTCA catheter through the gauze to remove any excess contrast medium.

4. Tighten the rotating hemostatic valve.
5. Reconfirm angiographic result.

NOTE: If an adequate expansion of the vessel has not been obtained, repeat the PTCA procedure (go back to point 9.3) or exchange to another balloon catheter or to a stent system of appropriate diameter.

6. Using a technique of choice, remove the PTCA

catheter, guidewire and guiding catheter from the vasculature.

10 IN VITRO INFORMATION

The data of **balloon compliance** charts printed on the product labels are based on *in vitro* testing. The **rated burst pressure** for the balloon dilatation is printed on the package label for the specific model. *In vitro* testing has shown that with 95% confidence, 99.9% of the balloons will not burst at or below the rated burst pressure. Balloons should not be inflated beyond the rated burst pressure.

11 REFERENCES

The physician should consult recent literature and current medical practice on balloon dilatation such as published by ACC/AHA (American Congress of Cardiology / American Heart Association).

12 DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY

This product was designed, manufactured, inspected, packaged, and sterilized with all reasonable care. The conditions under which the product is used, diagnosis and treatment of the patient, surgical procedure and / or its handling and storage after the product has left our possession directly affect the device and the results obtained from its use. The labelling and this instructions for use are expressly to be considered as an integral part of this disclaimer.

Acrostak's obligation under this warranty is limited to the replacement of the device but Acrostak shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. This warranty replaces all other warranties, whether expressed or implied, written or oral.

These limitations of liability and warranty are not intended to contravene any mandatory provisions of law applicable in the respective country. If any clause of the disclaimer is considered by a competent court to be invalid or to be in conflict with the applicable law, the remaining part of it shall not be affected and remain in full force and effect. The invalid clause shall be substituted by a valid clause which best reflects Acrostak's interest in limiting its liability or warranty.

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