## Cardiac Microcatheter

# M·CATH FLEXy

|                           |                  |    |
|---------------------------|------------------|----|
| INSTRUCTIONS FOR USE      | English          | 2  |
| GEBRAUCHSANWEISUNG        | Deutsch          | 4  |
| MODE D'EMPLOI             | Français         | 6  |
| ISTRUZIONI PER L'USO      | Italiano         | 8  |
| INSTRUCCIONES DE USO      | Español          | 10 |
| INSTRUÇÕES DE UTILIZAÇÃO  | Português        | 12 |
| GEbruIKSAANWIJZINGEN      | Nederlands       | 14 |
| BRUGSANVISNING            | Dansk            | 16 |
| BRUKSANVISNING            | Norsk            | 18 |
| BRUKSANVISNING            | Svenska          | 20 |
| INSTRUKCJA UŻYWANIA       | Polski           | 22 |
| NÁVOD K POUŽITÍ           | Česky            | 24 |
| NÁVOD NA POUŽITIE         | Slovenčina       | 26 |
| HASZNÁLATI UTASÍTÁS       | Magyar           | 28 |
| INSTRUCȚIUNI DE UTILIZARE | Română           | 30 |
| ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ            | Ελληνικά         | 32 |
| РУКОВОДСТВО ПО ПРИМЕНЕНИЮ | РУССКИЙ          | 34 |
| PETUNJUK PEMAČAIAN        | Bahasa Indonesia | 36 |
| GRAPHICAL SYMBOLS         |                  | 38 |

**Acrostak (Schweiz) AG**  
Stegackerstrasse 14  
8409 Winterthur, Switzerland  
Phone +41 52 233 95 51, [clinical@acrostak.com](mailto:clinical@acrostak.com)  
[www.acrostak.com](http://www.acrostak.com)

**QAdvis EAR AB**  
Ideon Science Park  
Scheelevägen 17  
SE-223 70 Lund, Sweden

**CE 2797**

## 1 DEVICE NAME

The device brand name is **M-CATH Flexy**. The generic device name is cardiac microcatheter.

## 2 DEVICE DESCRIPTION

**M-CATH Flexy** consists of a flexible polymer shaft with a stainless-steel braiding, a low-profile tapered tip, and a hub. The outer surface of the catheter is coated with a hydrophilic polymer having a high lubricity upon moistening, hence, enhancing manoeuvrability and patient comfort.

The microcatheter has one distal radiopaque marker near the distal tip. The shaft has a central guidewire lumen which permits the use of a guidewire with maximum diameter of 0.014" (0.36 mm). The guidewire lumen begins at the proximal luer port and ends at the distal tip.

**M-CATH Flexy** is introduced over a 0.014" guidewire with a 5F/0.056" guiding catheter.

|                               |                   |
|-------------------------------|-------------------|
| Maximum guidewire OD          | 0.36 mm (0.014")  |
| Minimum guiding catheter size | 5 F/0.056"        |
| Maximum pressure              | 2070 kPa (300psi) |
| Tensile strength              | 5 N               |

## 3 INTENDED PURPOSE

The **M-CATH Flexy** microcatheter is intended for supporting the guidewire in crossing stenotic lesions while performing PCI (percutaneous coronary intervention) and to facilitate the exchange of guidewires while maintaining access to distal vasculature.

## 4 MEDICAL CONDITIONS / INDICATIONS

The **M-CATH Flexy** microcatheter is indicated for:

- ☐ Coronary artery stenosis
- ☐ Bypass graft stenosis
- ☐ ST-elevation myocardial infarction (STEMI)

## 5 PATIENT POPULATION

The **M-CATH Flexy** microcatheter should only be used in adult persons above 18 years of age requiring percutaneous coronary intervention.

## 6 TARGET USER GROUP

**M-CATH Flexy** microcatheter should only be used by physicians, familiar in coronary angiography and PCI technique.

## 7 CONTRAINDICATIONS

- ☐ Patients with previous coronary artery spasm
- ☐ Not to be used for high pressure injections or in neurovasculature
- ☐ Not to be used for delivery of embolic materials.
- ☐ Not to be used in institutions where emergency coronary artery bypass graft surgery (CABG) cannot be performed.
- ☐ Not to be used in patients for whom the procedure is considered unsuitable by the physician.

## 8 WARNINGS

- ☐ This product is intended for single use only. Do NOT reprocess, resterilise and/or reuse. Processing and/or resterilisation measures can harm the characteristics of the catheter due to chemical, thermal and mechanical stresses. Furthermore, after reprocessing, there is an increased risk of contamination of the product as well as infection or cross-infection of the patient.
- ☐ Take care while handling and torquing the catheter during a procedure to reduce the possibility of accidental breakage, separation or kinking of the product.

- ☐ During the procedure, the catheter should only be manipulated under fluoroscopy.
- ☐ Never advance or withdraw the microcatheter against resistance until the cause of the resistance is determined by fluoroscopy. If not, the microcatheter and/or the guidewire could be damaged.
- ☐ Do not turn the microcatheter in the same direction, either clockwise or counterclockwise, for more than 5 consecutive turns. While the catheter tip is adjacent to the stenotic lesion. If resistance is observed, stop the movement and determine the cause of the resistance under fluoroscopy. Otherwise, the catheter or the vessels could be damaged, endangering the life of the patient.
- ☐ Do not turn the microcatheter excessively if it is bent, as this could break or damage the product.
- ☐ Do not use the catheter after the "Use by" date specified on the package.

## 9 PRECAUTIONS

- ☐ Ensure that the catheter is suitable for the intended procedure.
- ☐ The catheter should only be used by physicians trained in the performance of PCI.
- ☐ Utilise appropriate anticoagulant/antiplatelet and coronary vasodilator therapy during the procedure.
- ☐ During withdrawal of the microcatheter, hold a saline-soaked gauze around the exposed catheter shaft and pull the catheter through the gauze to remove any excess contrast medium.
- ☐ Prior to insertion or withdrawal of the microcatheter, wipe the guidewire with saline-soaked gauze to remove any contrast medium.
- ☐ Do not bend or kink the shaft or the attached catheter's hub.
- ☐ Do not expose the catheter to organic solvents (e.g., alcohol) or oleaginous contrast media. This could damage the product or cause loss of lubricity.
- ☐ There is no data available on the use of **M-CATH Flexy** during pregnancy or lactation. For safety reasons, pregnant women and breastfeeding mothers should therefore not be treated with **M-CATH Flexy**.

## 10 ADVERSE EFFECTS AND RISKS OF USE

These can include, but are not limited to:

- ☐ Dissection or perforation of the coronary artery
- ☐ Injury, rupture, intimal tear or other damage to the coronary artery
- ☐ Thrombosis, sudden and total occlusion of the coronary artery or bypass graft
- ☐ Unstable angina
- ☐ Arterial spasm
- ☐ Hypo/hypertension
- ☐ Hemorrhage or hematoma
- ☐ Allergic reactions
- ☐ Infections
- ☐ Acute myocardial infarction
- ☐ Arrhythmia
- ☐ Ventricular fibrillation
- ☐ Embolism
- ☐ Myocardial ischemia
- ☐ Death

## 11 CONTENTS

STERILE: Sterilized with ethylene oxide gas. Non-pyrogenic. Do not use if sterile package is opened or damaged.

### CONTENTS:

- One (1) Microcatheter
- One (1) Instructions for use manual

STORAGE: Store at room temperature, in a dry, dark place.

DISPOSAL INSTRUCTIONS: After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

## 12 DIRECTIONS FOR USE

### 12.1 Inspection Prior to Use

Carefully inspect the sterile package before opening. Do not use if any defects or kinks are noted or after the "Use by" date. If the integrity of the sterile packaging has been compromised before exceeding the product's "Use by" date, contact your local representative for return information.

### 12.2 Materials required (not included in device package)

- ☐ Appropriate guiding catheter(s) (see product label for dimensions)
- ☐ Syringe for flushing
- ☐ Heparinized normal saline solution (sterile)  $\leq 0.014$  inch (0.36 mm) guidewire

### 12.3 Preparation

1. Carefully remove the product within its dispenser from its packaging pouch.
2. Fill the dispenser slowly with heparinised saline solution using a syringe to soak the product. Check that the dispenser has been correctly filled with the solution and that the product is thoroughly wet.
3. Carefully remove the catheter from its dispenser.
4. Prime the product with the heparinised saline solution through the catheter hub in order to remove the air from inside.

### 12.4 Procedure

1. Prepare the vascular access site according to standard PCI practice.
2. Attach a haemostasis valve to the luer port of the guiding catheter positioned within the vasculature.
3. Insert the guidewire into the guiding catheter through the haemostasis valve.

**NOT** Make sure that the haemostasis valve is open

**E:** enough, otherwise it may cause resistance while inserting the product which could damage the device.

4. Under fluoroscopy, position the guidewire across the stenosis according to accepted PCI techniques.
5. Insert the proximal end of the guidewire into the distal tip of the microcatheter.

**NOT** Take care when inserting the microcatheter

**E:** into the haemostasis valve in order to avoid kinking/breakage.

6. Advance the microcatheter to the distal end of the guiding catheter.
7. Continue under fluoroscopy and use the radiopaque marker band to position the device close to the stenotic lesion.
8. To ensure that the guidewire has crossed the stenosis, a radiopaque dye should be injected. Make sure that the guidewire is correctly positioned using angiograms from multiple angles.
9. Keep the guidewire and the guiding catheter firmly in place, and slowly advance the distal end of the product gradually over the guidewire until the radiopaque marker has reached the desired position.

**NOT** To avoid product damage do not turn the catheter if it seems stuck or blocked.

10. The user can turn the microcatheter when necessary, during the intervention, but excessive turning of the microcatheter should be avoided. If resistance is observed, stop the movement, and determine the cause of the resistance as it could damage the catheter and/or the vessel.
11. Take care when re-inserting an angled guidewire into the product which remains in the coronary artery. Just before the tip of the guidewire protrudes from the tip of the product, stop advancing the guidewire and carefully withdraw the product so that the tip of the guidewire gradually protrudes from the product under high resolution fluoroscopy. The tip of the guidewire may spring out of the tip of the product, which may result in damage to the vessels.

#### 12.5 Catheter removal procedure

1. Fully open haemostasis valve.
2. While maintaining the guidewire position, withdraw the catheter.

**NOT** If you feel resistance, remove the catheter, guiding catheter and guidewire together.

3. Tighten the rotating haemostasis valve.
4. Reconfirm angiographic result.
5. Using a selected technique, remove the catheter, guidewire and guiding catheter from the vasculature.

**NOTE:**Wipe the product surface with heparinised saline solution and flush the inner lumen in order to remove residual blood.

**NOTE:**Handle the product carefully so as not to damage or kink the catheter.

#### 13 WARRANTY

Acrostak (Schweiz) AG guarantees that this product has been manufactured, packaged, and sterilised in compliance with Good Manufacturing Practice (GMP) for medical devices. Each product is individually tested before it is released for packaging.

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