



Short track to cross

ACROSS CTO™ **ST**

RX

- High pushability
- Strong 1mm tip
- Good lesion entry
- Low crossing profile
- Distal marker




Tip marker

Short 1 mm tip
Seamless Tip-to-shoulder

High RBP 17 Atm

Low entry
profile 0.016"



Ordering Information / Technical Information

Reference	Balloon Diameter	Balloon Length	Marker	Tip Length	Entry Profile	Max. Guidewire	Usable Length	Min. Guide Catheter / ID
110050483	1.1 mm	5 mm	1	1 mm	0.016"	0.014"	138 cm	5 F / 0.056"
110100483	1.1 mm	10 mm	1	1 mm	0.016"	0.014"	138 cm	5 F / 0.056"
110150483	1.1 mm	15 mm	1	1 mm	0.016"	0.014"	138 cm	5 F / 0.056"
110200483	1.1 mm	20 mm	1	1 mm	0.016"	0.014"	138 cm	5 F / 0.056"

Control your Size

Balloon Ø	Atm	Atm*	Atm	Atm	Atm	Atm	Atm**	Atm	Atm
	5	7	9	11	13	15	17	19	21
1.1 mm	1.04 mm	1.10 mm	1.13 mm	1.16 mm	1.20 mm	1.23 mm	1.28 mm	1.38 mm	1.45 mm

Atm*

Nominal Pressure (NP)

Atm**

Rated Burst Pressure (RBP)



The balloon material is a semi-compliant polyamide, providing a 'Work Your Size'.

Short track to cross