



High Performance in reaching and crossing

ACROSS HP™

RX

- Very low compliance
- Tapered tip for superior crossability
- Low crossing profile
- Outstanding material strength




Double marker

Seamless Tip-to-shoulder

High RBP up to 20 Atm

Low entry
profile 0.017"



Ordering Information / Technical Information

Reference	Balloon Diameter	Balloon Length	Marker	Tip Length	Entry Profile	Max. Guidewire	Usable Length	Min. Guide Catheter / ID
150100350	1.5 mm	10 mm	2	2 mm	0.017"	0.014"	138 cm	5 F / 0.056"
150150350	1.5 mm	15 mm	2	4 mm	0.017"	0.014"	138 cm	5 F / 0.056"
150200350	1.5 mm	20 mm	2	4 mm	0.017"	0.014"	138 cm	5 F / 0.056"
200100350	2.0 mm	10 mm	2	2 mm	0.017"	0.014"	138 cm	5 F / 0.056"
200150350	2.0 mm	15 mm	2	4 mm	0.017"	0.014"	138 cm	5 F / 0.056"
200200350	2.0 mm	20 mm	2	4 mm	0.017"	0.014"	138 cm	5 F / 0.056"
250100350	2.5 mm	10 mm	2	2 mm	0.017"	0.014"	138 cm	5 F / 0.056"
250150350	2.5 mm	15 mm	2	4 mm	0.017"	0.014"	138 cm	5 F / 0.056"
250200350	2.5 mm	20 mm	2	4 mm	0.017"	0.014"	138 cm	5 F / 0.056"
300100350	3.0 mm	10 mm	2	2 mm	0.017"	0.014"	138 cm	5 F / 0.056"
300150350	3.0 mm	15 mm	2	4 mm	0.017"	0.014"	138 cm	5 F / 0.056"
300200350	3.0 mm	20 mm	2	4 mm	0.017"	0.014"	138 cm	5 F / 0.056"
350100350	3.5 mm	10 mm	2	2 mm	0.017"	0.014"	138 cm	5 F / 0.056"
350150350	3.5 mm	15 mm	2	2 mm	0.017"	0.014"	138 cm	5 F / 0.056"
350200350	3.5 mm	20 mm	2	2 mm	0.017"	0.014"	138 cm	5 F / 0.056"
400100350	4.0 mm	10 mm	2	2 mm	0.017"	0.014"	138 cm	5 F / 0.056"
400150350	4.0 mm	15 mm	2	2 mm	0.017"	0.014"	138 cm	5 F / 0.056"
400200350	4.0 mm	20 mm	2	2 mm	0.017"	0.014"	138 cm	5 F / 0.056"
450100350	4.5 mm	10 mm	2	2 mm	0.017"	0.014"	138 cm	5 F / 0.056"
450150350	4.5 mm	15 mm	2	2 mm	0.017"	0.014"	138 cm	5 F / 0.056"
450200350	4.5 mm	20 mm	2	2 mm	0.017"	0.014"	138 cm	5 F / 0.056"

Control your Size



The balloon material is a semi-compliant polyamide, providing a 'Work Your Size'.

Balloon Ø	Atm 5	Atm 7	Atm 9	Atm* 11	Atm 13	Atm 15	Atm** 17	Atm** 18	Atm** 19	Atm** 20	Atm 21	Atm 22
1.5 mm	1.40 mm	1.43 mm	1.47 mm	1.50 mm	1.55 mm	1.60 mm	1.65 mm	1.67 mm	1.69 mm	1.71 mm	1.73 mm	1.75 mm
2.0 mm	1.82 mm	1.88 mm	1.94 mm	2.00 mm	2.07 mm	2.13 mm	2.18 mm	2.20 mm	2.22 mm	2.24 mm	2.25 mm	2.29 mm
2.5 mm	2.28 mm	2.34 mm	2.43 mm	2.50 mm	2.57 mm	2.65 mm	2.71 mm	2.74 mm	2.76 mm	2.78 mm	2.79 mm	2.81 mm
3.0 mm	2.65 mm	2.76 mm	2.88 mm	3.00 mm	3.10 mm	3.18 mm	3.24 mm	3.27 mm	3.30 mm	3.33 mm	3.35 mm	
3.5 mm	3.19 mm	3.31 mm	3.44 mm	3.50 mm	3.56 mm	3.62 mm	3.69 mm	3.74 mm	3.79 mm	3.83 mm		
4.0 mm	3.63 mm	3.71 mm	3.85 mm	4.00 mm	4.09 mm	4.18 mm	4.24 mm	4.28 mm	4.31 mm	4.34 mm		
4.5 mm	3.98 mm	4.15 mm	4.32 mm	4.50 mm	4.66 mm	4.77 mm	4.87 mm	4.91 mm	4.95 mm			

Atm*

Nominal Pressure (NP)

Atm**

Rated Burst Pressure (RBP)

High Performance in reaching and crossing